

ACCURACY YOU CAN ACT ON

Onyx® Vantage 9590

Fingertip Pulse Oximeter



Blue, black and red
color options



Carrying Case



Fingertip Lanyard



Fingertip Retractable
Lanyard

When you need to make decisions quickly, you do not have time to question the accuracy of your finger pulse oximeter. You can trust Nonin Medical's Onyx® Vantage finger pulse oximeter with PureSAT® technology to quickly and accurately capture SpO₂ and pulse rate measurements – even on patients where low perfusion is a challenge.¹ Accuracy matters, and Onyx Vantage provides accurate data for actionable decision making.¹

Reliable Readings for Confident Decisions¹



Accurate measurements across patient populations, inclusive of all skin pigmentations^{2,3}



Nonin PureSAT[®] technology captures SpO₂ and pulse rate measurements quickly and accurately¹



Withstands a minimum of 50 drops and exceeds IP32 water ingress testing requirements⁴



Clinician-facing, bright LED allows SpO₂ and pulse rate readings from any angle, in all light conditions



Fits pediatric to adult patients; accommodates fingers, thumbs or toes from 8mm to 25.4 mm (0.3 to 1.0 in.)



Three-color pulse quality indicator LED provides quick assessment of patient pulse quality

Designed for Lasting Value

Reliable readings in even the most challenging conditions¹

Power turns on when applied and off when removed

Small and portable, weighs less than two ounces

Provides up to 6,000 spot checks on two AAA batteries

Includes an industry-leading 4-year warranty

Onyx[®] Vantage
9590 Fingertip



Pulse quality indicator LED

Automatic On/Off

1. Nonin Medical, Inc. Clinical Data on File.

2. Bickler, P. MD, et. Al. (2005) *Effects of Skin Pigmentation on Pulse Oximeter Accuracy at Low Saturation*. Anesthesiology Vol. 102, pp 715-9

3. Bickler PE, Feiner JR, Severinghaus JW (2007). *Dark Skin Decreases the Accuracy of Pulse Oximeters at Low Oxygen Saturation: The Effect of Oximeter Probe Type and Gender*, Anesthesiology, December;105(6): S18-S23.

4. <https://www.youtube.com/watch?v=OYEXFYISJys>

Nonin.com

United States
info@nonin.com
+1-763-553-9968

International
infointl@nonin.com
+31 (0)13-45 87 130



FOCUS ON WHAT MATTERS MOST



WristOx₂® Model 3150 with Bluetooth® Low Energy
Wrist-Worn Pulse Oximeter



Wireless Connectivity



Compact Design



Long Battery Life

Featuring Nonin PureSAT® SpO₂ technology, the WristOx₂ 3150 with Bluetooth® Low Energy (BLE) delivers fast, accurate oxygen saturation and pulse rate readings in challenging conditions, including low perfusion and motion¹. The WristOx₂ 3150 with BLE is compatible with a wide range of Nonin PureLight® sensors and provides continuous, reliable SpO₂ monitoring, and data recording². Applications include oxygen qualification and sleep diagnostic screenings conducted at home, in the hospital, and in sleep labs.

For Proven Performance, Count on Nonin.



Memory Volume Indicator (MVI) Mode provides a graphic display, verifying that a sufficient amount of data has been recorded



Data can be transferred via USB connection to nVISION® or wirelessly to the NoninConnect app



1,080 hours of non-volatile memory at four-second sampling rate



Dependable design and construction, backed by a 3-year warranty



Exclusive sensors pair easily with the WristOx₂ 3150 oximeter for optimal patient comfort



Long lasting, 44-hour battery life operating on two AAA batteries

Experience Nonin Quality and Greater Affordability

Small size wears like a watch with large display screen and enhanced sensor connector design

Transfer data securely and wirelessly to the NoninConnect app*, available on the Apple Store or download via USB cable

Single unit—includes sensor of choice and wristband starter kit—includes sensor of choice, three wristbands, carry case, USB download cable, and nVISION® software

WristOx₂ Model 3150 with Bluetooth® Low Energy



WristOx₂ Accessories



8000S-WO2 Reusable Sensor

*The NoninConnect app isn't cleared as a medical device

1,2. Nonin Medical, Inc. Data on File.

Nonin.com

United States
info@nonin.com
+1-763-553-9968

International
info@nonin.com
+31 (0)13-45 87 130

